

llowing Products for Science™ www.biotium.com

Revised: May 24, 2024

Product Information

DAPI

| Cat. No. | Product | Unit Size | MW |
|----------|---|-----------|-----|
| 40009 | DAPI, dilactate | 10 mg | 457 |
| 40011 | DAPI, dihydrochloride | 10 mg | 350 |
| 40043 | DAPI, dilactate in H ₂ O, 10 mg/mL | 1 mL | 457 |

Storage and Handling

Store DAPI (solid form) desiccated at 4°C, protected from light. Store DAPI in $\rm H_2O$ at 4°C, protected from light. Product is stable for at least one year from date of receipt when stored as recommended.

Molecular Information

C₂₂H₂₇N₅O₆ (DAPI, dilactate)

C₁₆H₁₇Cl₂N₅ (DAPI, dihydrochloride)

CAS number

DAPI, dilactate: 28718-91-4 DAPI, dihydrochloride: 28718-90-3

Color and form: Yellow solid (40009, 40011); yellow liquid (40043)

Solubility: Soluble in water

Absorption/Emission: 352/464 nm (with DNA)

2CI-

$$H_2N_C$$
 H_2N_1
 H_2N_1
 H_2N_1

Figure 1. DAPI (4',6-Diamidino-2-Phenylindole), dihydrochloride

$$^{+}H_{2}N$$
 $^{+}H_{2}N$
 $^{+}H_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{+}NH_{2}$
 $^{-}NH_{2}$
 $^{-}NH_{2}$
 $^{-}NH_{2}$
 $^{-}NH_{2}$
 $^{-}NH_{2}$
 $^{-}NH_{2}$
 $^{-}NH_{2}$
 $^{-}NH_{2}$

Figure 2. DAPI (4',6-Diamidino-2-Phenylindole), dilactate

Product Description

DAPI is a blue DNA dye that is widely used as a nuclear counterstain for fluorescence microscopy, chromosome staining, and flow cytometry. The dye binds to the minor groove of dsDNA with approximately 20-fold fluorescence enhancement, with higher affinity for A-T rich regions. Because of the potential toxicity of DAPI, we offer DAPI as a ready-to-use solution in water as a safer alternative to weighing out the solid form. DAPI dilactate solid is more water soluble than the dihydrochloride salt of the dye.

At typical staining concentrations (~1 ug/mL) DAPI is impermeant to live cells, but useful as a nuclear counterstain in fixed cells or tissue sections. For live cell staining, we recommend using Hoechst cell-permeant blue DNA dyes (see Related Products).

DAPI and Hoechst have been shown to undergo UV-induced photoconversion, which causes fluorescence in both the blue channel and the FITC or Cy®3 channels. This could interfere with the imaging of other fluorescent probes or labeled antibodies. Biotium's NucSpot® Nuclear Stains (see Related Products) were designed as a convenient alternative to DAPI or Hoechst by offering nuclear-specific staining in 9 different colors, from green to near-IR, for flexible multiplexing.

References

1) Biochemistry 31, 3103 (1992); 2) Biochem Biophys Res Commun 170, 270 (1990); 3) J. Histochem Cytochem 38, 1323 (1990); 4) Methods Enzyme 168, 741 (1989); 5) Biochemistry 26, 4545 (1987); Biotechnic Histochem 70, 220 (1995).




Figure 3. Normalized excitation and emission spectra of DAPI with DNA.

Staining Protocols

Staining of fixed cells or tissue sections

- Dilute DAPI to 1 ug/mL in PBS. DAPI can be included together with antibodies or other probes, and can be diluted in buffers with detergent or blocking agents if convenient.
- 2 Add the staining solution to cells or tissue sections and incubate at room temperature for at least 5 minutes.
- 3 Image the samples with UV (355 nm) excitation in the DAPI channel.

Notes:

- a. Samples can be stored at 4°C after staining and before imaging.
- b. DAPI can be included directly in antifade mounting medium for one-step mounting and staining (see Biotium's EverBrite Mounting Media with DAPI under Related Products). When using DAPI in mounting medium, longer incubation times may be required for DAPI to completely penetrate the cell nuclei.

Staining bacteria or yeast

In bacteria, DAPI staining is dim compared to what is observed in mammalian cells. Live or killed bacteria can be stained with 10 ug/mL DAPI in PBS or 150 mM NaCI for 30 minutes at room temperature. DAPI tends to stain dead cells more brightly than live cells.

In *S. cerevisiae*, DAPI preferentially stains dead yeast with nuclear and cytoplasmic staining when used at 10 ug/mL in PBS; in live yeast DAPI shows dim mitochondrial staining.

Related Products

| Cat. No. | Product | | |
|-----------------|--|--|--|
| 40083 41040 | NucSpot® Nuclear Stains | | |
| 40044 | Hoechst 33258, 10 mg/mL in H ₂ O | | |
| 40045 | Hoechst 33258, pentahydrate | | |
| 40046 | Hoechst 33342, 10 mg/mL in H ₂ O | | |
| 40047 | Hoechst 33342, trihydrochloride trihydrate | | |
| 40081 | NucSpot® Live 488 | | |
| 40082 | NucSpot® Live 650 | | |
| 40085 | NucSpot® Far-Red | | |
| 40060 | RedDot™1 Far-Red Nuclear Stain | | |
| 40061 | RedDot™2 Far-Red Nuclear Stain | | |
| 40084 | 7-AAD Solution, 1 mg/mL | | |
| 40048 | Propidium Iodide in Buffer, 50 ug/mL | | |
| 23001 | EverBrite™ Mounting Medium | | |
| 23002 | EverBrite™ Mounting Medium with DAPI | | |
| 23003 | EverBrite™ Hardset Mounting Medium | | |
| 23004 | EverBrite™ Hardset Mounting Medium with DAPI | | |
| 23008 | Drop-n-Stain EverBrite™ Mounting Medium | | |
| 23009 | Drop-n-Stain EverBrite™ Mounting Medium with DAPI | | |
| 23016 | EverBrite™ Hardset Mounting Medium with NucSpot® 640 | | |
| 30068 | ViaFluor® 405 SE Cell Proliferation Kit | | |
| 30086 | ViaFluor® 488 SE Cell Proliferation Kit | | |
| 70065 | LipidSpot™ 488 Lipid Droplet Stain | | |
| 70069 | LipidSpot™ 610 Lipid Droplet Stain | | |
| 70066 | LysoView™ 405 Lysosome Stain | | |
| 70067 | LysoView™ 488 Lysosome Stain | | |
| 70061 | LysoView™ 540 Lysosome Stain | | |
| 70058 | LysoView™ 633 Lysosome Stain | | |
| 70059 | LysoView™ 650 Lysosome Stain | | |
| 70070 | MitoView™ 405 Mitochondrial Stain | | |
| 70054 | MitoView™ Green Mitochondrial Stain | | |
| 70055 | MitoView™ 633 Mitochondrial Stain | | |
| 70075 | MitoView™ 650 Mitochondrial Stain | | |
| 70068 | MitoView™ 720 Mitochondrial Stain | | |
| 70064 | ViaFluor® 405 Live Cell Microtubule Stain | | |
| 70062 | ViaFluor® 488 Live Cell Microtubule Stain | | |
| 70063 | ViaFluor® 647 Live Cell Microtubule Stain | | |
| 30090 | CellBrite® Fix 488 Membrane Stain | | |
| 30088 | CellBrite® Fix 555 Membrane Stain | | |
| 30089 | CellBrite® Fix 640 Membrane Stain | | |
| 30092- 30099 | MemBrite® Fix Cell Surface Staining Kits | | |
| 22023 | Paraformaldehyde, 4% in PBS, Ready-to-Use Fixative | | |

Please visit our website at www.biotium.com for information on our life science research products, including novel live cell stains for nuclei, mitochondria, and other organelles, fluorescent CF® Dye antibody conjugates and reactive dyes, apoptosis reagents, fluorescent probes, and kits for cell biology research.

Cy Dye is a registered trademark of Cytiva.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.